

Flacons de pesée en verre borosilicaté 3.3



- couvercle rodé
- conformes DIN 12605

capacité	dim. (Ø x h)	rodage	référence	Prix HT
rodage 19/12				
3 ml	Ø 25 x 25 mm	19/12	TK1025	
rodage 24/12				
30 ml	Ø 30 x 70 mm	24/12	TK1033	
rodage 34/12				
45 ml	Ø 40 x 60 mm	34/12	TK1043	
66 ml	Ø 40 x 80 mm	34/12	TK1045	
rodage 45/12				
20 ml	Ø 40 x 30 mm	45/12	TK1046	
54 ml	Ø 50 x 50 mm	45/12	TK1048	
88 ml	Ø 50 x 70 mm	45/12	TK1049	
rodage 50/12				
126 ml	Ø 60 x 70 mm	50/12	TK1055	
150 ml	Ø 60 x 80 mm	50/12	TK1056	



- couvercle rodé
- conformes DIN 12605

capacité	dim. (Ø x h)	rodage	référence	Prix HT
rodage 19/12				
3 ml	Ø 25 x 25 mm	19/12	TK2025	
5 ml	Ø 25 x 30 mm	19/12	TK2026	
9 ml	Ø 25 x 40 mm	19/12	TK2027	
13 ml	Ø 25 x 50 mm	19/12	TK2028	
rodage 24/12				
7 ml	Ø 30 x 30 mm	24/12	TK2029	
13 ml	Ø 30 x 40 mm	24/12	TK2030	
19 ml	Ø 30 x 50 mm	24/12	TK2031	
24 ml	Ø 30 x 60 mm	24/12	TK2032	
30 ml	Ø 30 x 70 mm	24/12	TK2033	
36 ml	Ø 30 x 80 mm	24/12	TK2034	
rodage 29/12				
10 ml	Ø 35 x 30 mm	29/12	TK2035	
18 ml	Ø 35 x 40 mm	29/12	TK2036	
26 ml	Ø 35 x 50 mm	29/12	TK2037	
34 ml	Ø 35 x 60 mm	29/12	TK2038	
42 ml	Ø 35 x 70 mm	29/12	TK2039	
rodage 34/12				
13 ml	Ø 40 x 30 mm	34/12	TK2040	
23 ml	Ø 40 x 40 mm	34/12	TK2041	
34 ml	Ø 40 x 50 mm	34/12	TK2042	
45 ml	Ø 40 x 60 mm	34/12	TK2043	
55 ml	Ø 40 x 70 mm	34/12	TK2044	
66 ml	Ø 40 x 80 mm	34/12	TK2045	
rodage 45/12				
20 ml	Ø 50 x 30 mm	45/12	TK2046	
37 ml	Ø 50 x 40 mm	45/12	TK2047	
54 ml	Ø 50 x 50 mm	45/12	TK2048	
88 ml	Ø 50 x 70 mm	45/12	TK2049	
105 ml	Ø 50 x 80 mm	45/12	TK2050	
rodage 50/12				
29 ml	Ø 60 x 30 mm	50/12	TK2051	
53 ml	Ø 60 x 40 mm	50/12	TK2052	
78 ml	Ø 60 x 50 mm	50/12	TK2053	
102 ml	Ø 60 x 60 mm	50/12	TK2054	
126 ml	Ø 60 x 70 mm	50/12	TK2055	
150 ml	Ø 60 x 80 mm	50/12	TK2056	

Flacons de pesée en verre



- couvercle rodé
- conformes DIN 12605

capacité	dim. (Ø x h)	référence	Prix HT
Pèse-filtres, forme basse 1			
5 ml	Ø 28 x 25 mm	GW1280	
15 ml	Ø 38 x 30 mm	GW1281	
30 ml	Ø 54 x 30 mm	GW1282	
80 ml	Ø 85 x 30 mm	GW1283	
Pèse-filtres, forme haute			
10 ml	Ø 28 x 40 mm	GW1285	
20 ml	Ø 32 x 50 mm	GW1286	
45 ml	Ø 38 x 70 mm	GW1287	
70 ml	Ø 44 x 80 mm	GW1288	